

SKILL[®]

Surgical Instruments





SKILL® The World's Best Surgical Instruments

All **SKILL®** Instruments are manufactured on a three-pillar principle, namely raw material quality, skilled craftsmanship as well as precise and validated heat treatment and tempering.

Our goal is achieved by using German Forgings and German Steel for the production of our **SKILL®** Instruments. Using only German Craftsmanship, engineering and know-how throughout the manufacturing process along with the validated heat treatment is a necessity in the production of our **SKILL®** Instruments.

SKILL® Instruments are manufactured in Germany according to the EU Directive 93/42/EWG, as well as the following ISO Standards;

ISO: 9001
ISO:13458
ISO: 7153-1

No compromise is made on the quality of a **SKILL®** Surgical Instrument.
The **SKILL®** short chain of supply allows clients to purchase the highest quality of Instrument at the most competitive pricing.

SKILL® Instruments offer the Surgeon a heightened sense of Instrument control allowing the Surgeon to greatly reduce the trauma and risk associated with Surgery.

OPERATE WITH SKILL®, OPERATE WITH US.

SKILL® Instruments carry a Lifetime warranty which covers any material defect or faulty manufacturing.

SKILL® Instruments will be repaired or replaced free of charge without compromise to the product in the event of a valid claim.

SKILL® Instruments are to be returned to the authorised distributor for maintenance, any work carried out by an unauthorised person will void the warranty.






No part of the Catalogue may be copied, forwarded or distributed without the express written consent given by the Manufacturer or its Distributor.

Operate with **SKILL®**
Surgical Instruments

Index

	Page
Table	4
Preface	5
Information	6
1/1 Size Containers	7
1/1 Size Containers	Basic Model Standard Model 8
1/1 Size Containers	Safe Model 9
1/1 Size Containers	Bottoms and Lids 10
1/1 Size Containers	Wire Basket 11
3/4 Size Containers	12
3/4 Size Containers	Basic Model 13
3/4 Size Containers	Safe Model 14
3/4 Size Containers	Bottoms and Lids 15
3/4 Size Containers	Wire Basket 16
1/2 Size Containers	17
1/2 Size Containers	Basic Model Standard Model 18
1/2 Size Containers	Safe Model 19
1/2 Size Containers	Bottoms and Lids 20
1/2 Size Containers	Wire Basket 21
Other Containers	22
Containers and Bottoms	23
Lids and Accessories	24
Mini Containers	25
Bottoms, Lids and Accessories	26-27
Implant Containers	28
Scopy Containers	28
Containers and Bottoms	29
Removal Containers	Lids 30
Removal Containers	Stainless Steel and Aluminium 31
Accessories	32
General Containers	Filters, Indicators and Labels 33
General Containers	Silicon Mats and Fixation Clamps 34
Accessories	Components 35
Permanent Filters	36
Silicon Fixation	37



Sterilization container	Packing		Steam sterilization procedure acc. to DIN 58946, Part I				
	Description	Symbolic representation and packing system	Gravity procedure	Pre-vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injection procedure
With filters (STB - F)	Textiles in sterilization container with filter in lid and bottom acc. to DIN 58952 Part I		+	+	+	+	+
			0	0	+	+	+
	Textiles in sterilization container with filter in lid and bottom not perforated acc. to DIN 58952 Part I		—	—	+	+	+
			—	—	+	+	0
	Textiles in sterilization container with filter in bottom lid, not perforated acc. to DIN 58952 Part I		0	+	+	+	+
			—	0	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III in sterilization container with filter in lid and bottom acc. to DIN 58952 part I		+ ³⁾	+	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III in sterilization container with filter in lid and bottom not perforated acc. to DIN 58952 part I		—	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket acc. to DIN 58952 part III in sterilization container with filter in bottom acc. to DIN 58952 part I		0	+	+	+	+
With valves (STB - V)	Textiles in sterilization container with valves in lid and bottom acc. to DIN 58952 Part I		—	0	+	+	0
			—	—	+	+	0
	Textiles in sterilization container with valves in lid no valves on the bottom		—	—	+	+	—
			—	—	+	+	—
	Instruments in wire basket in sterilization container with valve in lid and bottom		—	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket in sterilization container with valve in lid and not perforated		—	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
Explanation of symbols: + applicable - not applicable 0 limited application							
1. In order to avoid mistakes with regard to service operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures. 2. Only sterilization containers guaranteeing vapo ration or discharge of condensate are suitable above. 3. Only sterilization containers guaranteeing vapo ration or discharge of condensate are suitable above.							

Preface

The **SKILL** sterilization container systems allow you a systematic organisation of your whole instrument sterilization.

The comfortable handling and the diversity of our container give you safety and control to your instrument sterilization. Stocking, transport and disposal are in a perfect organisation with our container.

Our sterilization container systems are according to the requirements of international standards.

SKILL sterilization container are made of anodised aluminium, have little weight and very good drying properties.

Colour coded lids and identification labels allow you the exact control of your instrument sterilization.

Our wide product range keeps you a proper adjustment to your needs.



Choose between our standard and our safety model.

The safety lid protects against contamination while transport and stocking.

Filter- and safety lid can only separate by opening container.





The locking flaps are movable 180°
on a pivot and can be secured in
any position.



The filter retainer is locked by pressing down in center of the filter.
Opening the filter retainer by pressing the plate


1/1 Container

1/1 Size Containers



1/1 Size Containers

Basic Model | Standard Model




	FF-060-01
	FF-060-02
	FF-060-03
	FF-060-04
	FF-060-05
	FF-060-06

Weight
2634 g



Outside	580 x 280 x 100 mm
Inside	550 x 265 x 84 mm



	FF-061-01
	FF-061-02
	FF-061-03
	FF-061-04
	FF-061-05
	FF-061-06

Weight
3156 g

	FF-060-11
	FF-060-12
	FF-060-13
	FF-060-14
	FF-060-15
	FF-060-16

Weight
3129,5 g



Outside	580 x 280 x 135 mm
Inside	550 x 265 x 115 mm

	FF-061-11
	FF-061-12
	FF-061-13
	FF-061-14
	FF-061-15
	FF-061-16

Weight
3244 g

	FF-060-21
	FF-060-22
	FF-060-23
	FF-060-24
	FF-060-25
	FF-060-26

Weight
3219 g



Outside	580 x 280 x 150 mm
Inside	550 x 265 x 130 mm

	FF-061-21
	FF-061-22
	FF-061-23
	FF-061-24
	FF-061-25
	FF-061-26

Weight
3382 g

	FF-060-31
	FF-060-32
	FF-060-33
	FF-060-34
	FF-060-35
	FF-060-36

Weight
3429,5 g



Outside	580 x 280 x 200 mm
Inside	550 x 265 x 180 mm

	FF-061-31
	FF-061-32
	FF-061-33
	FF-061-34
	FF-061-35
	FF-061-36

Weight
3738,5 g

	FF-060-41
	FF-060-42
	FF-060-43
	FF-060-44
	FF-060-45
	FF-060-46

Weight
3919 g



Outside	580 x 280 x 260 mm
Inside	550 x 265 x 240 mm


	FF-061-41
	FF-061-42
	FF-061-43
	FF-061-44
	FF-061-45
	FF-061-46

Weight
4272,5 g

1/1 Size Containers

Safe Model

Bottom non perforated




	FF-062-01
	FF-062-02
	FF-062-03
	FF-062-04
	FF-062-05
	FF-062-06

Weight
3486 g



Outside	580 x 280 x 100 mm
Inside	550 x 265 x 84 mm

Botto perforated



	FF-063-01
	FF-063-02
	FF-063-03
	FF-063-04
	FF-063-05
	FF-063-06

Weight
4007,5 g

	FF-062-11
	FF-062-12
	FF-062-13
	FF-062-14
	FF-062-15
	FF-062-16

Weight
3981,5 g



Outside	580 x 280 x 135 mm
Inside	550 x 265 x 115 mm

	FF-063-11
	FF-063-12
	FF-063-13
	FF-063-14
	FF-063-15
	FF-063-16

Weight
4095,5 g

	FF-062-21
	FF-062-22
	FF-062-23
	FF-062-24
	FF-062-25
	FF-062-26

Weight
4070,5 g



Outside	580 x 280 x 150 mm
Inside	550 x 265 x 130 mm

	FF-063-21
	FF-063-22
	FF-063-23
	FF-063-24
	FF-063-25
	FF-063-26

Weight
4233 g

	FF-062-31
	FF-062-32
	FF-062-33
	FF-062-34
	FF-062-35
	FF-062-36

Weight
4280 g



Outside	580 x 280 x 200 mm
Inside	550 x 265 x 180 mm

	FF-063-31
	FF-063-32
	FF-063-33
	FF-063-34
	FF-063-35
	FF-063-36

Weight
4589 g

	FF-062-41
	FF-062-42
	FF-062-43
	FF-062-44
	FF-062-45
	FF-062-46

Weight
4770,5 g





Outside	580 x 280 x 260 mm
Inside	550 x 265 x 240 mm

	FF-063-41
	FF-063-42
	FF-063-43
	FF-063-44
	FF-063-45
	FF-063-46

Weight
5140 g


1/1 Size Containers

Bottoms and Lids



	Outside	Inside	
FF-064-10	580 x 280 x 100 mm	550 x 262 x 84 mm	1923 g
FF-064-20	580 x 280 x 135 mm	550 x 262 x 117 mm	2019 g
FF-064-30	580 x 280 x 150 mm	550 x 262 x 132 mm	2108 g
FF-064-40	580 x 280 x 200 mm	550 x 262 x 182 mm	2317 g
FF-064-50	580 x 280 x 260 mm	550 x 262 x 244 mm	2808 g

	Outside	Inside	
FF-065-10	580 x 280 x 100 mm	550 x 262 x 84 mm	2044 g
FF-065-20	580 x 280 x 135 mm	550 x 262 x 117 mm	2132 g
FF-065-30	580 x 280 x 150 mm	550 x 262 x 132 mm	2270 g
FF-065-40	580 x 280 x 200 mm	550 x 262 x 182 mm	2627 g
FF-065-50	580 x 280 x 260 mm	550 x 262 x 244 mm	2851 g



580 x 280 mm

Grey	FF-066-01
Yellow	FF-066-02
Orange	FF-066-03
Red	FF-066-04
Blue	FF-066-05
Green	FF-066-06

Weight
1111 g




Grey	FF-067-01
Yellow	FF-067-02
Orange	FF-067-03
Red	FF-067-04
Blue	FF-067-05
Green	FF-067-06


Weight
851,5 g

1/1 Size Containers


Wire Basket




Wire Basket		
FF-069-01	485 x 255 x 30 mm	1108 g
FF-069-02	485 x 255 x 50 mm	1112 g
FF-069-03	485 x 255 x 70 mm	1176 g
FF-069-05	485 x 255 x 100 mm	1289 g



Wirebasket		
FF-069-11	540 x 255 x 30 mm	1149 g
FF-069-12	540 x 255 x 50 mm	1189 g
FF-069-13	540 x 255 x 70 mm	1259 g
FF-069-14	540 x 255 x 100 mm	1354 g



Wire Base		
FF-069-21	545 x 255 x 30 mm	729 g



Drape Retainer		
FF-965-12	545 x 255 mm	



Drape		
FF-965-07	130 x 60 cm	
FF-965-08	140 x 100 cm	
FF-965-09	185 x 45 cm	
FF-965-10	185 x 60 cm	


3/4 Size Containers



3/4 Size Containers

Basic Model

Bottom non perforated



	FF-160-01
	FF-160-02
	FF-160-03
	FF-160-04
	FF-160-05
	FF-160-06

Weight
2268,5 g



Outside	465 x 280 x 100 mm
Inside	425 x 260 x 84 mm

Bottom perforated



	FF-161-01
	FF-161-02
	FF-161-03
	FF-161-04
	FF-161-05
	FF-161-06

Weight
2503 g

	FF-160-11
	FF-160-12
	FF-160-13
	FF-160-14
	FF-160-15
	FF-160-16

Weight
2549,5 g



Outside	465 x 280 x 135 mm
Inside	425 x 260 x 115 mm

	FF-161-11
	FF-161-12
	FF-161-13
	FF-161-14
	FF-161-15
	FF-161-16

Weight
2715 g

	FF-160-21
	FF-160-22
	FF-160-23
	FF-160-24
	FF-160-25
	FF-160-26

Weight
2620,5 g



Outside	465 x 280 x 150 mm
Inside	425 x 260 x 130 mm


	FF-161-21
	FF-161-22
	FF-161-23
	FF-161-24
	FF-161-25
	FF-161-26

Weight
2718 g

3/4 Size Containers

Safe Model

Bottom non perforated




	FF-162-01
	FF-162-02
	FF-162-03
	FF-162-04
	FF-162-05
	FF-162-06

Weight
2980 g



Outside	465 x 280 x 100 mm
Inside	425 x 260 x 84 mm

Bottom perforated



	FF-163-01
	FF-163-02
	FF-163-03
	FF-163-04
	FF-163-05
	FF-163-06

Weight
3213,5 g

	FF-162-11
	FF-162-12
	FF-162-13
	FF-162-14
	FF-162-15
	FF-162-16

Weight
3260 g



Outside	465 x 280 x 135 mm
Inside	425 x 260 x 115 mm

	FF-163-11
	FF-163-12
	FF-163-13
	FF-163-14
	FF-163-15
	FF-163-16

Weight
2425 g

	FF-162-21
	FF-162-22
	FF-162-23
	FF-162-24
	FF-162-25
	FF-162-26

Weight
2975 g



Outside	465 x 280 x 150 mm
Inside	425 x 260 x 130 mm


	FF-163-21
	FF-163-22
	FF-163-23
	FF-163-24
	FF-163-25
	FF-163-26

Weight
3561 g


3/4 Size Containers

Bottoms and Lids

Bottom non perforated



Bottom perforated



	Outside	Inside	
FF-164-01	465 x 280 x 100 mm	425 x 262 x 85 mm	1396 g
FF-164-02	465 x 280 x 135 mm	425 x 262 x 121 mm	1676 g
FF-164-03	465 x 280 x 150 mm	425 x 262 x 135 mm	1792 g

	Outside	Inside	
FF-165-01	465 x 280 x 100 mm	425 x 262 x 85 mm	1629 g
FF-165-02	465 x 280 x 135 mm	425 x 262 x 121 mm	1842 g
FF-165-03	465 x 280 x 150 mm	425 x 262 x 135 mm	1978 g

Lid perforated



465 x 280 mm

Safety Lid



	FF-166-01
	FF-166-02
	FF-166-03
	FF-166-04
	FF-166-05
	FF-166-06

Weight
873 g

	FF-167-01
	FF-167-02
	FF-167-03
	FF-167-04
	FF-167-05
	FF-167-06


Weight
710 g

3/4 Size Containers

Wire Basket



Wire Basket		
FF-169-01	405 x 255 x 30 mm	972 g
FF-169-02	405 x 255 x 50 mm	988 g
FF-169-03	405 x 255 x 70 mm	1061 g
FF-169-04	405 x 255 x 100 mm	1143 g



Wire Base		
FF-169-11	405 x 255 x 30 mm	598 g



Drape Retainer	
FF-965-13	405 x 255 mm




Drape	
FF-165-08	140 x 100 cm

1/2 Size Containers



1/2 Size Containers

Basic Model | Standard Model




	FF-260-01
	FF-260-02
	FF-260-03
	FF-260-04
	FF-260-05
	FF-260-06

Weight
1744,5 g



Outside	285 x 280 x 100 mm
Inside	262 x 260 x 84 mm



	FF-261-01
	FF-261-02
	FF-261-03
	FF-261-04
	FF-261-05
	FF-261-06

Weight
1890,5 g

	FF-260-12
	FF-260-13
	FF-260-14
	FF-260-15
	FF-260-16
	FF-260-17

Weight
1947 g



Outside	285 x 280 x 135 mm
Inside	262 x 260 x 116 mm

	FF-261-11
	FF-261-12
	FF-261-13
	FF-261-14
	FF-261-15
	FF-261-16

Weight
2125 g

	FF-260-21
	FF-260-22
	FF-260-23
	FF-260-24
	FF-260-25
	FF-260-26

Weight
1981 g



Outside	285 x 280 x 150 mm
Inside	262 x 260 x 132 mm

	FF-261-21
	FF-261-22
	FF-261-23
	FF-261-24
	FF-261-25
	FF-261-26

Weight
2462,5 g

	FF-260-31
	FF-260-32
	FF-260-33
	FF-260-34
	FF-260-35
	FF-260-36

Weight
2101 g



Outside	285 x 280 x 200 mm
Inside	262 x 260 x 182 mm

	FF-261-31
	FF-261-32
	FF-261-33
	FF-261-34
	FF-261-35
	FF-261-36

Weight
2513,5 g

	FF-260-41
	FF-260-42
	FF-260-43
	FF-260-44
	FF-260-45
	FF-260-46

Weight
2483,5 g



Outside	285 x 280 x 260 mm
Inside	262 x 260 x 238 mm


	FF-261-41
	FF-261-42
	FF-261-43
	FF-261-44
	FF-261-45
	FF-261-46

Weight
2657,5 g

1/2 Size Containers

Safe Model

Bottom non perforated




	FF-262-01
	FF-262-02
	FF-262-03
	FF-262-04
	FF-262-05
	FF-262-06

Weight
2162,5 g



Outside	285 x 280 x 100 mm
Inside	262 x 260 x 84 mm

Bottom perforated



	FF-263-01
	FF-263-02
	FF-263-03
	FF-263-04
	FF-263-05
	FF-263-06

Weight
2309,5 g

	FF-262-11
	FF-262-12
	FF-262-13
	FF-262-14
	FF-262-15
	FF-262-16

Weight
2351 g



Outside	285 x 280 x 135 mm
Inside	262 x 260 x 116 mm

	FF-263-11
	FF-263-12
	FF-263-13
	FF-263-14
	FF-263-15
	FF-263-16

Weight
2466 g

	FF-262-21
	FF-262-22
	FF-262-23
	FF-262-24
	FF-262-25
	FF-262-26

Weight
2394 g



Outside	285 x 280 x 150 mm
Inside	262 x 260 x 132 mm

	FF-262-21
	FF-262-22
	FF-262-23
	FF-262-24
	FF-262-25
	FF-262-26

Weight
2583 g

	FF-262-31
	FF-262-32
	FF-262-33
	FF-262-34
	FF-262-35
	FF-262-36

Weight
2560 g



Outside	285 x 280 x 200 mm
Inside	262 x 260 x 182 mm

	FF-263-31
	FF-263-32
	FF-263-33
	FF-263-34
	FF-263-35
	FF-263-36

Weight
2785,5 g

	FF-262-31
	FF-262-32
	FF-262-33
	FF-262-34
	FF-262-35
	FF-262-36

Weight
2930 g



Outside	285 x 280 x 260 mm
Inside	262 x 260 x 238 mm


	FF-263-41
	FF-263-42
	FF-263-43
	FF-263-44
	FF-263-45
	FF-263-46

Weight
3069,5 g


1/2 Size Containers

Bottoms and Lids

Bottom non perforated



Bottom perforated



	Outside	Inside	
FF-264-01	285 x 280 x 100 mm	262 x 262 x 82 mm	1096 g
FF-264-02	285 x 280 x 135 mm	262 x 262 x 120 mm	1284 g
FF-264-03	285 x 280 x 150 mm	262 x 262 x 135 mm	1327 g
FF-264-04	285 x 280 x 200 mm	262 x 262 x 185 mm	1447 g
FF-264-05	285 x 280 x 260 mm	262 x 262 x 240 mm	1817 g

	Outside	Inside	
FF-265-01	285 x 280 x 100 mm	262 x 262 x 82 mm	1242 g
FF-265-02	285 x 280 x 135 mm	262 x 262 x 120 mm	1398 g
FF-265-03	285 x 280 x 150 mm	262 x 262 x 135 mm	1471 g
FF-265-04	285 x 280 x 200 mm	262 x 262 x 185 mm	1672 g
FF-265-05	285 x 280 x 260 mm	262 x 262 x 240 mm	2000 g

Lid perforated



285 x 280 mm

Safety Lid



	FF-262-01
	FF-262-02
	FF-262-03
	FF-262-04
	FF-262-05
	FF-262-06


Weight
654 g

	FF-267-01
	FF-267-02
	FF-267-03
	FF-267-04
	FF-267-05
	FF-267-06


Weight
458 g

1/2 Size Containers


Wire Baskets



Wire Basket		
FF-269-01	255 x 245 x 30 mm	637,5 g
FF-269-02	255 x 245 x 50 mm	648,5 g
FF-269-03	255 x 245 x 70 mm	691,5 g
FF-269-04	255 x 245 x 100 mm	758,5 g



Wire Base		
FF-269-11	258 x 258 x 30 mm	351 g



Drape Retainer	
FF-965-17	255 x 255 mm



Drape	
FF-265-04	90 x 60 cm
FF-265-05	90 x 80 cm
FF-265-06	130 x 45 cm
FF-265-07	130 x 60 cm


Other Containers

Wire Baskets



Flat Containers

Containers and Bottoms




	FF-360-01
	FF-360-02
	FF-360-03
	FF-360-04
	FF-360-05
	FF-360-06


Weight
1285 g




Outside	285 x 280 x 55 mm
Inside	270 x 270 x 33 mm




	Outside	Inside	
FF-362-01	285 x 280 x 55 mm	270 x 270 x 33 mm	611 g
FF-362-02	285 x 280 x 85 mm	270 x 270 x 64 mm	703 g



	Outside	Inside	
FF-362-03	285 x 280 x 85 mm	270 x 270 x 64 mm	857 g



	Outside	Inside	
FF-362-01	285 x 280 x 55 mm	270 x 270 x 33 mm	611 g
FF-362-02	285 x 280 x 85 mm	270 x 270 x 64 mm	703 g



	Outside	Inside	
-362-03	285 x 280 x 85 mm	270 x 270 x 64 mm	857 g

Flat Containers

Lids and Accessories




	FF-362-11
	FF-362-12
	FF-362-13
	FF-362-14
	FF-362-15
	FF-362-10


Weight
674 g




Drape Retainer	
FF-965-17	255 x 255 mm



Wire Basket	
FF-269-01	255 x 245 x 30 mm
FF-269-02	255 x 245 x 50 mm



Wire Base		
FF-269-11	245 x 255 x 30 mm	351 g




Silicone Mat	
FF-965-21	245 x 245 mm



Drape	
FF-965-04	90 x 60 cm
FF-965-05	90 x 80 cm

Flat Containers

Mini Containers




	FF-460-01
	FF-460-02
	FF-460-03
	FF-460-04
	FF-460-05
	FF-460-06

Weight
510,5 g



Outside	300 x 140 x 40 mm
Inside	280 x 130 x 27 mm



	FF-461-01
	FF-461-02
	FF-461-03
	FF-461-04
	FF-461-05
	FF-461-06

Weight
589 g

	FF-460-11
	FF-460-12
	FF-460-13
	FF-460-14
	FF-460-15
	FF-460-16

Weight
567 g



Outside	300 x 140 x 70 mm
Inside	280 x 130 x 55 mm

	FF-461-11
	FF-461-12
	FF-461-13
	FF-461-14
	FF-461-15
	FF-461-16

Weight
630 g

	FF-460-21
	FF-460-22
	FF-460-23
	FF-460-24
	FF-460-25
	FF-460-26

Weight
625,5 g




Outside	300 x 140 x 100 mm
Inside	280 x 130 x 84 mm

	FF-461-21
	FF-461-22
	FF-461-23
	FF-461-24
	FF-461-25
	FF-461-26


Weight
625 g

Mini Containers


Bottoms, Lids and Accessories



	Outside	Inside	
FF-462-01	300 x 140 x 40 mm	280 x 130 x 27 mm	190 g
FF-462-02	300 x 140 x 70 mm	280 x 130 x 55 mm	241 g
FF-462-03	300 x 140 x 100 mm	280 x 130 x 84 mm	305 g



	Outside	Inside	
FF-462-11	300 x 140 x 40 mm	280 x 130 x 27 mm	258 g
FF-462-12	300 x 140 x 70 mm	280 x 130 x 55 mm	304 g
FF-462-13	300 x 140 x 100 mm	280 x 130 x 84 mm	326 g



	FF-462-21
	FF-462-22
	FF-462-23
	FF-462-24
	FF-462-25
	FF-462-26

Weight
325,5 g



Drape Retainer	
FF-965-18	280 x 125 mm




Silicone Mat	
FF-965-20	270 x 125 mm



Drape	
FF-965-01	60 x 40 cm

Dental Containers

Bottoms, Lids and Accessories



FF-5602-01

Weight
783 g



310 x 190 x 40 mm
280 x 182 x 35 mm

FF-560-02

Weight
862 g




310 x 190 x 65 mm
280 x 182 x 55 mm

FF-560-03


Weight
1022 g



310 x 190 x 130 mm
280 x 182 x 124 mm



	Outside	Inside	
FF-560-11	310 x 190 x 40 mm	280 x 182 x 35 mm	341 g
FF-560-12	310 x 190 x 65 mm	280 x 182 x 55 mm	415 g
FF-560-13	310 x 190 x 130 mm	280 x 182 x 124 mm	570 g



	Lid	
FF-560-21	310 x 190 mm	417 g



	Drape Retainer
FF-965-19	310 x 190mm



	Drape
FF-965-01	60 x 40 cm
FF-965-02	65 x 65 cm
FF-965-03	85 x 85 cm

Implant Containers



Lid perforated · Bottom perforated (non Filter)



	FF-660-11
	FF-660-12
	FF-660-13
	FF-660-14
	FF-660-15
	FF-660-16

Outside	500 x 155 x 75 mm
Inside	495 x 155 x 75 mm

Weight
654 g




	FF-660-01
	FF-660-02
	FF-660-03
	FF-660-04
	FF-660-05
	FF-660-06

Outside	500 x 155 x 75 mm
Inside	495 x 155 x 75 mm

Weight
368 g

Scopy Containers



Lid perforated · Bottom perforated (non Filter)



	FF-660-22
--	-----------

Outside	450 x 70 x 70 mm
Inside	450 x 70 x 60 mm

Weight
302 g



	FF-660-21
--	-----------

Outside	250 x 60 x 50 mm
Inside	250 x 55 x 50 mm


Weight
152 g

Removal Containers

Containers and Bottoms



	FF-760-01
	FF-760-02
	FF-760-03
	FF-760-04
	FF-760-05
	FF-760-06



Outside	285 x 280 x 100 mm
Inside	262 x 260 x 82 mm

Weight
1677 g



	FF-760-11
	FF-760-12
	FF-760-13
	FF-760-14
	FF-760-15
	FF-760-16

Outside	285 x 280 x 150 mm
Inside	262 x 260 x 132 mm

Weight
1929 g




	FF-760-21
	FF-760-22
	FF-760-23
	FF-760-24
	FF-760-25
	FF-760-26

Outside	580 x 280 x 150 mm
Inside	548 x 258 x 130 mm

Weight
2883 g



Bottom non perforated



FF-064-30	580 x 280 x 150 mm	2108 g
-----------	--------------------	--------

FF-264-01	285 x 280 x 150 mm	1096 g
FF-264-03	285 x 280 x 150 mm	1327 g

Removal Containers

Lids



285 x 280 mm

	FF-760-51
	FF-760-52
	FF-760-53
	FF-760-54
	FF-760-55
	FF-760-56

Weight
458 g



580 x 280 mm

	FF-760-41
	FF-760-42
	FF-760-43
	FF-760-44
	FF-760-45
	FF-760-46

Weight
875 g

Removal Containers

Stainless Steel and Aluminium



	HE-939-20	For 20 brushes
	HE-939-10	For 10 brushes



Brush container
Stainless steel
Brushes are taken out by rotating the handle under the container.
Container are supplied together with wall mounting fixtures.

	FF-860-01
	FF-860-02
	FF-860-03
	FF-860-04
	FF-860-05

Brush container. For 20 brushes
240 x 275 x 70 mm
Aluminium
Brushes are taken out by rotating the handle under the container.
Container are supplied together with wall mounting fixtures.



Accessories



General Containers

Filters, Indicators and Labels



FF-960-01
Paper Filter
VE 100 St.
Container
1/1, 3/4 und 1/2
Pack of 100
For 1/1, 3/4 and
Half Size Containers



FF-960-02
Textilfilter
Container
1/1, 3/4 und 1/2
Single textile
For 1/1, 3/4 and Half
Size Containers

FF-960-03
Disposable Paper Filter
Pack of 100
For Mini Containers
FF-960-05
100 Pcs/Package
for Dental Containers

FF-960-04
Textilfilter / Single textil
For Mini Containers
FF-960-06
Textilfilter / Single textil
for Dental Containers



FF-960-10
Pack of 100
Paper label with indicator
For 1/1, 3/4 and Half Size Containers



FF-960-11
Pack of 100
Paper label with indicator
For Mini and Dental Containers



Filter tray FF-960-08

For Paper Filter
FF-960-09 For Textile
Filter




FF-960-07
Filter Dispenser



FF-960-12 Plastic
Security Seal Pack
of 100



FF-960-29
Dispenser for identification labels



	FF-960-21
	FF-960-22
	FF-960-23
	FF-960-24
	FF-960-25
	FF-960-26


Coloured identification labels.
Printed to your specification, max.16 letters

General Containers


Silicon Mats and Fixation Clamps



FF-965-20	270 x 125 mm
-----------	--------------



FF-965-21	245 x 245 mm
-----------	--------------



FF-965-24	220 x 230 mm
FF-965-25	380 x 230 mm
FF-965-23	440 x 230 mm
FF-965-26	520 x 230 mm



FF-960-50
24 x 28 x 25
ø 7 - 12 mm



FF-960-51
24 x 28 x 29
ø 12 - 16 mm



FF-960-52
24 x 28 x 31
ø 16 - 26 mm



FF-960-53
24 x 28 x 40
ø 26 - 36 mm



FF-960-54
24 x 28 x 44
ø 36 - 45 mm



FF-960-48
50 x 20 x 12 mm



FF-690-49
65 x 20 x 20 mm



Accessories


Components



FF-960-31
Holding Pin for Partition Sheets



FF-960-32
Holding Pin for Wire Baskets



Sheets	
FF-960-33	50 x 20 mm
FF-960-34	130 x 20 mm
FF-960-35	225 x 20 mm
FF-960-36	460 x 20 mm



Sliding Rails with Securing Screws
ER841.000 230 mm
ER841.010 470 mm



FF-960-44 230 mm
FF-960-45 470 mm

Limiting Bar with Securing Screws



FF-960-41
Divider, to be used
with FF-960-39



FF-960-46
Spiral Holding




FF-960-47
Double Spiral Holding

Permanent Filters




FF-960-13

Life for 1200 Sterilization cycles
for 1/1, 3/4 and 1/2 Container



FF-960-14




Life for 1200 Sterilization cycles
for Dental Container



FF-960-15



Life for 1200 Sterilization cycles
for Mini Container



Silicone Fixation



Silicone fixation

The silicone fixation protects your instruments during sterilization, transport and storage. The division of the instruments in this type of fixation is flexible. We adapt the silicone fixation of instruments to your needs.





Authorised Dealer / Agent

www.skillinstruments.co.za

SKILL® Surgical Instruments
Rottweiler Street 53
D-78669 Wellendingen
Germany